Dartguns Remote Drug Delivery (RDDD) In Clinical Practice
Feedlot Management School
Saskatoon Sask Feb 12th

By Roy Lewis DVM



Roy Lewis Bio


- Roy Lewis DVM 1981 WCVM
- 1981-2012 Westlock Veterinary Center
- Focused on beef herd health and reproduction
- Large Animal practice 5%-10 bison & any species "chuteable"
- Owned up to 70 Bison cows
- Worked for the pharma industry Merck Animal Health Jan 2012-May 2020
- Part of Family Owned Purebred Cattle (Simmental & Angus) Grain & Potatoes (Lewis Farms)
- Always a good day when working with cattle or bison

Veterinarians Involvement & Responsibility

- Dartguns (RDD) are here to stay need to use wisely.
- Maximizing Safety (training) in usage and dart placement.
- Help with Animal Welfare for bison/cattle and producer
- Must use wisely following BQA VBP+ Guidelines
- Avoid Residues (Both meat from antimicrobials and metal)
- VBP+ Producers Must Record ID and Product Used.
- Collar may stay in. (foreign body) Biodegradable collar



When is Remote Drug Delivery Appropriate?




Scan for more RDDD resources, including a 4-minute video.



Find more practical, science-based information for beef cattle producers at BeefResearch.ca.

Remote Drug Delivery Devices (RDDD)







98% sold as ½ needles to go S.C. Sask Distributor 2019 30,000 packs 1 pack=5 darts(150,000 darts)







Safety First for Operator and Intended Target Practise Makes perfect Crossbows large volume darts not desirable.

Pneu Dart Scope vs Sights Holographic Sights






Pros of Usage


- Convenience (Less labour)
- Darting much easier skill than roping.
- Earliest treatment (treated from tractor, quad horse or feed truck!!
- Safety of both cattle & producer/veterinarian (ie treating bulls)
- Less Stress (animal welfare) Not roping and choking down
- Most time gets in right location



- 98% of sales ½ inch to get s.c.
- ¾ inch cows and bulls .
- Never seen a broken dart with pneu-dart
- Odd report of dart in the US
 embedded in the musculature
 These make headlines!
- Can see broken needles with the Dan-inject system. Another Dart system.
- Single use dart with Pneu Dart
- Use with tranquilizers

2 Canadian Distributors 1 Sask & 1 Alberta





Canadian Pneu Dart Distributors Veterinary Clinics & Farm Supply Stores 50/50 Split

AB

Central Alberta Ag Supply is the Canadian Distributor for Pneu-Dart. Please see the map to find Pneu-Dart Dealers across

For questions about Pneu-Dart please call our office and we will be happy to help you.

1 | Dougan's Farm & Feed Supply

403-638-3959 1015 Main Ave E Sundress's, AB ToM 1Xo 2 | Battlefords Animal Hospital

306-446-8387 11204 Railway Ave E. #7A North Battleford, SK S9A 2R7

3 | East Anglie Breeding Center Ltd.

403-362-3296 728-2 Ave E. Brooks, AB T1R 1B5 4 | Hackman Feeds

306786-7713 385 5th Ave N Yorktown, SK S3n2Z9

5 | Keddies Tack & Western Wear

780-532-4888 9808-132 Ave Grand Prairie, AB T8V 4J6

6 | Park Range Veterinary Services

306-764-6998 380A S Industrial Dr Prince Albert, SK S6V 7Ys

7 | Peace River Building Products

250-785-2116 9511-85th Ave Fort St John, BC V1J 6M5 8 | Sunny South Vet Clinic 403-327-2425

1118-41St N Lethbridge, AB T1H 5M1

9 | Total Farm Supply Ltd.

204-727-3276 1530-18 St N Brandon, MB R7C 1A6 10 | Weir Veterinary Services

MB

780-875-2281 5603-44St Lloydminster, AB T9V oB2

11 | Westlock Vet Centre

780-349-3663 10004-94a Ave Westlock, AB T7P 2M9 12 | Winfield Feed & Farm Supply

780-682-2493 Main St Winfield, AB ToC 2Xo 13 | Hooves & Paws Vet

306-278-2252 116 Main St Elrose, SK SoL oZo 14 | Barr North Vet Services


780-674-2100 6302 49 St Barrhead, AB T7N 1A3

Screenshot




Accuracy of Darting

- 5/15 missed their mark
- All still S.C. if ½ inch needle
- Scope increases accuracy
- Loss of some darts in the environment
- 90% of guns now sold with scopes .
 Holographic sites so easy to use.



Articles on Pros and Cons

- Doc Talk Dr. Dan Thomson Doc Talks 422
- Some training videos web sites.
- www.pneudart.com
- in USA 70% of cattle producers use a remote delivery system
- Location & Residues.
- BCRC in @2023 Survey, Brochure
- QR code leading to the 4 minute video
- Video Available on BCRC Website along with many other factual stuff.
- Dr John Campbell Podcast on Dart Usage



Mitigating Risks With RDDDs

- Try and Shoot from 10-30 ft
- Shoot from as close to 90 degrees as you can
- Make a decent distance exam first.
- Consider everything you use Extra label!!!

- Maintain your gun
- Watch expiry date of darts.
- Waxy collar dries out and falls off so there is an expiry date

Cons of Dart Usage


- Right Diagnosis (Footrot) Distance
 Exam
- Product not given by right route IM vs SC
- Muscle damage and inappropriate location. Volume
- Withdrawal times altered!!
- Safety muscle damage hit scapula,
- Could go IV (never been reported to date)
- Cost of equipment.

- Often using high category
 antimicrobials because low
 dosage, long acting (doesn't
 follow prudent usage guidelines.
- Beef Quality Assurance
- Darts left in muscle (can happen if charge too high very very rare)
- May not get ID of animal records VBP+ in field or feedlot??
- Use binoculars and go on









Lacerated Prepuce Good prognosis 70-80% with Time.



Underdosing NOT discharging fully

- 4/15 Coetze et al 2018
- 3/15 Rivera et al 2019
- Using products absolutely should not be given this way
- Micotil (safety for both operator and cattle (given IV by mistake) Elanco doesn't want it given in a RDDD
- Excede given not even close to the label. Not to be given SC or IM period other than base of ear.

Cons of Usage # 2

- drug residues
- Metal residues darts left in muscle (very very rare)



Tips to Darting

- Make sure clinical diagnosis is correct and it is truly necessary to dart.
- Don't over- power the dart.
- Check to make sure fully discharged when retrieve.
- Use neck and try and avoid butt muscles
- ½ inch dart all that is necessary except big bulls.
- Record on medical records which ones have been treated.



